HOOFCARE SPRAYMAT

HoofCare SprayMat for prevention and treatment of contagious hoof diseases

The hoofCare SprayMat will be installed on the floor of a walkway just outside the exit of the milking parlour or milking robot. If the cow walks over the mat and steps with one foot on or near a colored sphere, a valve will open and the liquid is then forced out through the nozzle.

No obstacle for the animals

The liquid is sprayed force fully against the foot and the medication or disinfecting agent will be sprayed into the hard-to-reach areas. The system is connected to the water mains. A dosing pump in the control unit will add the correct amount of medication to the water.

The Vink HoofCare SprayMat hits all the right spots!

- Medication is consistently applied to the correct areas
- Reduces labor costs
- Medication does not get Contaminated/always new clean liquid
- · No obstacle for the animals
- No electricity or electronics
- Provides an effective and economical use of consumables

The dimension of one mat is 80x40 cm, Behind a milking robot you will need 4 mats to reach a length of 160 cm. Behind a milking parlor we recommend 5 mats for a total length of 200 cm

The mat, which requires virtually no maintenance, can be left in a permanent position as it is not an obstacle to animals. The HoofCare SprayMat is ideally positioned at the exit of a parlour or milking robot where animals can cross it a number of times per day. Ensures all hoofs are disinfected several times per week with always clean liquid.





Location:

The hoof care spray mat consists of three special designed parts, bottom mat, spray mat and cover mat. The mat can be left in a permanent position as it is not an obstacle to animals, which requires no maintenance. The animals will cross the mat a number of times per day. Which ensures all hoofs are disinfected several times per week.







